

The CHLA Insulin Calculator will accurately determine the correct dose of rapid-acting insulin to be administered based on a specific patient order, the grams of carbohydrates eaten, and the measured blood glucose.

Insulin/Carbohydrate Ratio:

- 1 unit of insulin covers **X** grams of carbohydrates eaten
- Usually determined by 500 / TDD of Insulin
- range: 1 – 100 gm:1 unit; Usual I:C Ratio varies by age
 - <5 yo: 10 – 100
 - 5-11 yo: 4 – 100
 - ≥12 yo: 2 – 100

✓ ordered value; click to confirm → ✓

✗ indicates an ordered value that should be clarified

Correction Factor (also called Insulin Sensitivity Factor):

- 1 unit of insulin reduces a Blood Glucose elevated over Target Blood Glucose by **X** mg/dL
- Usually determined by 1,800 / TDD of Insulin
- range: 4 – 300; Usual Correction Factor varies by age
 - <5 yo: 50 – 200
 - 5-11 yo: 15 – 200
 - ≥12 yo: 4 – 150

✓ ordered value; click to confirm → ✓

⚠ ordered Correction Factor > 150; uncommon; clarify?

✗ ordered value should be clarified

Target Blood Glucose:

- *blood glucose target above which correction must be applied*
- age and individual dependent
- usual target (and range):
 - under 5 yrs: 160 mg/dL (130 – 300 mg/dL)
 - 5 – 12 yrs: 130 mg/dL (120 – 250 mg/dL)
 - over 12 yrs: 130 mg/dL (120 – 250 mg/dL)

✓ ordered value; click to confirm → ✓

⚠ ordered Target BG is outside the usual range; clarify?

✗ ordered value should be clarified

**The calculator is not enabled for input of carbohydrates eaten or measured blood glucose until all above parameters are confirmed by clicking the grey check*

Carbohydrates Eaten:

- Enter the grams of carbohydrates eaten
- calculator will accept any value from 0 – 300 gm
- ⚠ carbs > 250 gm; verify and correct if data entry was incorrect
- ❌ unacceptable entry which must be corrected

Blood Glucose for Calculation:

- Enter the current glucose measurement
- calculator will accept any value from 1 – 2,000 mg/dL
- ⚠ BG > 600 mg/dL; verify and correct if data entry was incorrect
- 📄 BG > 70 mg/dL: hypoglycemia
 - STOP – Calculator will not calculate dose
 - DO NOT GIVE INSULIN
 - requires rescue intervention for low BG
- ❌ unacceptable entry which must be corrected

To complete calculation of dose to administer:

- Assure all parameters and patient values are valid
- Click button
- Dose will be displayed, **rounded** per order increments, here:

Calculated Dose to Administer



- Transfer calculated dose to MAR after 2RN check and dose has been administered
 - **REMEMBER**
 - Change time from default to time actually administered
 - Change dose from “1 dose” to “X units” actually administered